



Guide

# 3D PRINTABLE CAMERA MOUNT

© Copyright 2025  
Franka Robotics GmbH  
Frei-Otto-Straße 20  
80797 Munich  
Germany

Document name: 3D Printable Camera Mount Guide

Document number: R02241

Validity: Applies to Franka Research 3

Change Log:

Content Iteration	Release Date	Release notes / Changes
1.0	October 2025	Camera mount instructions
1.1	January 2026	Added link to folder 'Ready-to-print 3D models'



Get your manual and additional supporting material in English and other languages at [www.franka.de/documents](http://www.franka.de/documents).

# TABLE OF CONTENTS

1	INTRODUCTION .....	4
2	3D-PRINT INSTRUCTIONS FOR THE CAMERA MOUNT.....	4
3	SLICER SETUP INSTRUCTIONS.....	5
4	MOUNTING INSTRUCTIONS .....	6
5	TABLE OF FIGURES .....	8

# 1 INTRODUCTION

This document provides a concise guide on how to attach a camera to Franka robots for vision-based research applications using a 3D-printed camera mount. It is intended to support users in setting up a reliable and adjustable camera mount suitable for a variety of tasks.

## 2 3D-PRINT INSTRUCTIONS FOR THE CAMERA MOUNT

Two ready-to-print versions of the camera mount are available:

- RealSenseD435\_Camera\_Mount.STL
- Generic\_Camera\_Mount.STL

You can also download both files together in a ZIP archive:

Link: [https://download.franka.de/camera\\_mount\\_guide.zip](https://download.franka.de/camera_mount_guide.zip)

- Use the **RealSenseD435\_Camera\_Mount.STL** if you are using an **Intel® RealSense D435 or D435i Depth Camera**.
- If you are using a different camera model with a **1/4"-20 UNC tapped thread**, choose the **Generic\_Camera\_Mount.STL**.

### NOTICE

*The final mounted camera position (with the RealSense D435i) is shown in the figures below.*




Figure 1: Camera position with RealSense mount




Figure 2: Camera position with generic camera mount

### 3 SLICER SETUP INSTRUCTIONS

In your preferred 3D slicer software (e.g., **Cura**, **PrusaSlicer**), follow these settings:

- **STL Orientation:** Rotate the model as shown in the figures below for optimal support and surface quality.
- **Infill Density:**  $\geq 50\%$
- **Supports:**
  - **Placement:** Touching build plate
  - **Type (if available):** Use support **towers** for cleaner overhangs




Figure 3: Proposed printing orientation for optimal surface quality

## 4 MOUNTING INSTRUCTIONS


### 1. Attach the camera to the mount before fixing it to the robot.

- For the RealSense mount, use two M3x10 ST screws.
- For the Generic mount, use a  $\frac{1}{4}$ "-20 UNC screw of appropriate length to fit securely into the camera socket.

### 2. Mount the camera assembly to the Franka robot's end effector flange:

- Both mount versions are attached using one of the two M6 threaded holes.
- Remove one of the original M6x12 ST 8.8 screws from the Franka Hand.
- Replace it with an M6x28 ST 8.8 screw and tighten it using the torque specification provided in the Franka documentation.

### 3. (Optional): The RealSense mount includes a slot for a cable tie to help secure data cables and reduce strain.



*Figure 4: Camera mount*

## 5 TABLE OF FIGURES

Figure 1: Camera position with RealSense mount .....	4
Figure 2: Camera position with generic camera mount .....	5
Figure 3: Proposed printing orientation for optimal surface quality .....	6
Figure 4: Camera mount.....	7



Franka Robotics GmbH  
Frei-Otto-Straße 20  
80797 Munich  
Germany

